

| PRODUCT | QUANTITY |
|---------|----------|
|---------|----------|

|                  |        |
|------------------|--------|
| TICO® Yellow 622 | 98,03% |
|------------------|--------|

|                  |       |
|------------------|-------|
| TICO® Orange 640 | 1,97% |
|------------------|-------|

Pigment level calculated on solid binder: 23,7%

Paint system: alkyd melamine resin based on Alkydal F300/Maprenal MF650 and other stoving enamel systems

| PRODUCT | QUANTITY |
|---------|----------|
|---------|----------|

|                  |        |
|------------------|--------|
| TICO® Yellow 622 | 98,70% |
|------------------|--------|

|                  |       |
|------------------|-------|
| TICO® Orange 640 | 1,30% |
|------------------|-------|

Pigment level calculated on solid binder: 26,9%

Paint system: 2K-Polyurethane, solvent based on Desmophen A365/Desmodur N3390